DESTINAZIONE D'USO

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www.nkf.it

INTENDED USE

For in vitro diagnostic use only. Reagents to be used on NIHON KOHDEN cell counters in combination as indicated in **Table 1**. For the intended use see **Table 2**. For professional use only.

STORAGE AND STABILITY AFTER FIRST OPENING

Store between 2° and 8 °C. It can be transported at 1 ° -30 ° C within 4 weeks. Install each reagent after checking the integrity of the packaging and the validity of the device. Once opened, the reagent may be used until depleted and by the date indicated in **Table 3**.

COMPOSITION AND TARGET POPULATION

Active ingredients See **Table 4**. Target population see **Table 5**.

WARNINGS

Labeling information See **Table 6**. Any serious incidents that may occur involving the device must be reported both to the manufacturer and to the competent authority of the Member State in which the incident has occurred.

Information on disposal: Product remnants, waste from its use and empty containers must be disposed of in accordance with applicable legislation. Waste from the instrument, such as potentially infected substances, must be handled with care and disposed of in accordance with the applicable legislation.

SAMPLE AND MARKER

The reagents were created for use with whole blood samples. Refer to the Operating Manual of the instrument for sample collection and storage procedures. Marker: reticulocyte in whole blood.

PROCEDURE AND MEASUREMENT PRINCIPLES

Please refer to the Operating Manual of the instrument for information on how to install the reagents, how to preflush the lines, and how to use the instrument. Operate at the operating temperatures indicated in the Operating Manual of the instrument. Use each product in its original packaging (do not decant). Operate according to good laboratory practice (GLP).

Refer to **Table 7** for the measurement principles.

INTERFERING SUBSTANCES OR LIMITATIONS

The following interfering substances have been confirmed to have no effect on reticulocyte measurement below the listed concentrations:

Substance	Max. serum concentration
Bilirubin C	21.1 mg/dL
Bilirubin F	21.1 mg/dL
Hemolysis hemoglobin	0.51 g/dL
Chyle	1490 FTU
Total protein	3.99 g/dL

BIBLIOGRAPHY

Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition, CLSI Document H26-A2, Vol. 30, No. 14, 2010.

METHODOLOGY AND CALCULATION OF THE ANALYTICAL RESULT

Refer to the Operating Manual of the instrument.

INTERNALI QUALITY CONTROL

Background check

After the instrument lines have been preflushed, check the background (background count) as explained in the Operating Manual of the instrument, and make sure that the values obtained for each parameter are within the limits specified in the Operating Manual. If one or more of the parameters are not within these limits, proceed according to the instrument manufacturer's instructions. If the problem can be attributed to one of the reagents, replace the package in use with an intact one and repeat the check. Do not proceed with the analysis until the background values are acceptable. Calibration and calibration verification procedures. It is recommended to check the calibration of the instrument under the following conditions:

1. after the initial installation of the instrument
2. at least every six months (or, in the case of more restrictive indications, as recommended by the manufacturer of the instrument)
3. after a technical intervention deemed critical (replacement of components that may affect instrument performance, software upgrade, complete replacement of reagents)
4. if quality control data reveals a significant change in baseline accuracy (drift), or the data is outside the established acceptability limits.

Before performing calibration procedures, it is recommended to carry out preventive maintenance (baseline accuracy and precision check: ref. CLSI document H26-A2). Calibration of the instrument and calibration checking must be performed using appropriate material (fresh whole blood samples with values that have been assigned using manual reference methods or an appropriate calibrator), following the instructions provided in the Operating Manual of the instrument. If the instrument needs to be calibrated, it is recommended that the calibration check be performed immediately after calibration, following the instructions in the Operating Manual of the instrument (also refer to CLSI document H26-A2).

Daily quality control procedures

In order to ensure the accuracy and reliability of the results obtained, it is recommended to perform the quality control procedures following the instructions in the Operating Manual of the instrument, using appropriate control material (appropriate stabilized control blood), immediately before, during and immediately after each analysis sequence. It is recommended to use two different levels of stabilized control blood (refer to CLSI document H26-A2).

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NIHON KOHDEN FIRENZE S.r.l.